



Developed By P-I Brånemark



MT-F

THE NEXT GENERATION



SMART GUIDE

MT-F

THE NEXT GENERATION



SIMPLE EXPERIENCE FOR EXCEPTIONAL OUTCOMES

Enhanced Biological Metrics to unlock immediate replacement potential

The P-I Implant Systems were developed by Professor Per-Ingvar Brånemark, the Osseointegration pioneer, jointly with scientists from renowned universities and the P-I Research & Development team to meet the modern implant dentistry demands.

In 2012, Ospol AB Sweden was acquired, and key technologies were integrated in the P-I solutions.

With the human biology, long-term expertise, clinical and scientific evidences as a foundation, our main objective is to support you in patient-focused treatments by providing Implant Systems that represent: Simplification • High Performance • Safety and Longevity

MT-F is the Next Generation System, a result of the P-I Brånemark fundamentals evolutionized by outstanding Biological Metrics and Simplicity.

Contents

MT-F UNIQUENESS	06
------------------------	----

P-I Technologies

 MT Interface	10
 OSPOL Surface	14
 ENHANCED Ceramic Surface	18
Custom Instruments	19


MT-F Implants	26
Soft Tissue Healing	29
Conical Abutment	31
Abutment Cemented Cylinder	33
Contour & Esthetic Abutments	35
Cylinder over Implant	37
 DIGITAL	39

Kit & Instruments	44
 Guided & Conventional SURGERY	47
Resonance Frequency	53
Locator, Overdenture	53
Surgical Sequence	54
Torques, Implant Packaging	57
LifeTime Guarantee	59



Adaptive bone contact

Multiple transitions • Interpolated core



Less bone displacement

In all bone densities

● MT-F

Less Trauma • Site Engagement

Cutting threads • Pronounced depth in all sections

● MT-F UNIQUENESS

The multiplicity of interpolated core transitions associated with the P-I Conical Drills site preparation, and the gradual evolution of the pronounced depth cutting threads, result in a gentle implant-to-osteotomy engagement in all sections independently.

These unique geometrical combinations provide greater initial contact area with significantly less bone displacement and compression enhancing the Biological Metrics.



↑ ISQ


ITV ↗

RTQ %

Less compression • Progressive torque • Greater area

MT-F displaces significantly less bone volume and achieves similar or higher Insertion Torque Value in all bone densities, exhibiting greater area in comparison to the leading competitive tapered-active implants of similar dimensions. Data on file.





○ MT-F

Enhanced Biological Metrics

↑ ISQ

ITV ↗

RTQ %

Biological Metrics

High initial and secondary Implant Stability Quotient • ISQ measurements by Resonance Frequency Analysis • RFA in association with sufficient Insertion Torque Value • ITV and low rotational micro-mobility, indicated by the proportional Removal Torque • RTQ% to the obtained ITV, are relevant Biological Metrics and critical success factors for the prosthetic rehabilitation of patients with implants in post extraction, healed and compromised sites, low density bone and in combination with tissue regeneration techniques.

The P-I expertise

Our expertise related to ISQ using RFA micromovement measurements to clinically monitor Osseointegration and to determine when to load implants, originates from the acquisition of Ospol AB in 2012.

Ospol AB and Osstell AB were sister companies established in Sweden and developers of an Implant System and RFA measurement technologies, respectively. The Ospol AB developments of the last 20 years are comprised within the P-I Implant Systems and the newest technologies are present in The Next Generation • MT-F System.



sister companies



Increased coronal space

Slightly inward flange

Cortical stability

Micro Patterns

Adaptive bone contact

Interpolated core transitions



Gradual thread evolution


Pronounced depth in all sections

Early engagement


Gentle cutting • Dual thread

Axial insertion control

Biological Width positioning



● MT-F



Peri-implant tissue preservation

● **MT Interface**

Biological
Width



Double
Sealing
System stability

Interface in clinical use for 15+ years • Superior biomechanics • Double Sealing

The P-I Morse Taper is an original technology. Highlighting 8.5° x 2 conical indexed, 3mm long, the P-I MT Interface offers a high torsional yield and fatigue strength when compared to other leading systems and was adopted by a global leader in 2015. The MT-F Ø 3.3 Implant can withstand static load of approximately 600N. Data on file.

The high-preload Double Sealing mechanism has easy prosthetic reversibility, seals the Abutment on the MT Interface and the MT Screw on the Abutment, stabilizing the system, minimizing micromovement and microleakage in comparison to certain leading systems under simulated occlusal stress. The Double Sealing is an important hypothesis for the clinical consideration of MT-F Implant placement observing Biological Width principles.

One Interface


Various prosthetic platforms



High Pre-Load

Effective sealing

P-I Coating



Easy reversibility

Low stress to peri-implant tissue



Sealing starts at provisionalization

P-I Coating+ is a biocompatible layer that reduces friction and, combined with the MT Screw optimized geometry, provides a substantially higher and homogeneous pre-load, clamping, in comparison to titanium screws at the same tightening torque of 25 Ncm. Data on file. MT Retriever is used to cancel the morse sealing and safely remove Abutments.



Prosthetic Overview

Screw-Retained

Cement-Retained

Screw • Cemented

D DIGITAL
prosthetics



N **R** **W**



D One MT Interface for all MT-F Implant dimensions



Locator® Abutments available



Same MT Screw and Prosthetic Driver Ø 1.2 for all Abutments, except straight Conical Abutment and Locator®. All P-I Components are supplied with the respective screw, except Milling Blank.

D See DIGITAL Prosthetics. Libraries are available at pibranemark.com and/or from CAD/CAM software system. Direct Geometries, 3D model Analogs and Geometry Contours available. Link C post has Dentsply Sirona, CEREC® dimensions. Always check latest library version and availability.



files for open
source software
pibranemark.com


3shape  exocad



Strong Osseointegration

REDUCTION OF BIOFILM INFECTIONS

BIOACTIVE



 **OSPOL**
SURFACE

Improved bone response

In comparison to rougher oxidized and blasted surfaces

Less bacterial adhesion

Equivalent to turned surfaces • Minimally rough

Chemically enhanced

Anodized • Bioactive ions

CP **Ti**
grade 4

Anodized
+
Chemistry


ZERO %
Aluminum
or Acid

Widely documented Evolution of moderately rough surfaces




Direct and early response

Oxide, micropores and crystal structures
greatly influence bone response



Biochemical bond, bone in-growth and mechanical interlocking



Courtesy of : YT Sul, A. Wennerberg,
T. Albrektsson

Surface chemistry, anodic oxidation and ion incorporation, enhance Osseointegration and compensate for minimal roughness

OSPOL Surface was developed in the Gothenburg University, Sweden, and is documented in several publications. In continual evolution since 2000 and in clinical use for over 15 years, the OSPOL Surface is a modern technology for a rapid and strong bone response and its modification method was adopted by a global leader in 2019. Less prone to bacterial adhesion, it is a pioneer technology for chemical modification of thin anodized, oxidized, ion incorporation of smoother implants surfaces.




Higher ISQ for chemically- modified Surface

OSPOL Surface modification method achieves faster secondary Implant Stability Quotient • ISQ measured by Resonance Frequency Analysis RFA indicating potential for shorter healing periods.

Less bacterial adhesion and biofilm formation

BIOACTIVE

⚠ OSPOL Surface is easier to clean than rougher surfaces and its bioactivity reduces biofilm formation. The bacterial adhesion is similar to turned, machined surfaces.



(!) Some conditions, whether combined or not, represent contraindications, limitations and risks, relative and absolute, for the treatment of patients with implants. There are several risk factors in Osseointegration widely described in literature. ISQ is a critical factor to clinically monitor Osseointegration. Data from pre-clinical studies.


ENHANCED Ceramic Surface



D DIGITAL
prosthetics

The P-I ENHANCED Ceramic Surface for Prosthetic Components is a smooth, tailored anodized surface technology that can potentially offer additional benefits for soft tissue in comparison to untreated titanium.

The acid-free P-I ENHANCED Ceramic Surface also exhibits a golden color for improved esthetics.



Kit


Surgical & Prosthetic



Stainless Steel

BIOSAFETY & DURABILITY
Conventional

Conical Drills



Less friction. Less trauma

Constant apical conical angle • 3 cutting areas

Corrosion protection


Wear resistance • Diamond Like Carbon

State-of-the-art

Special P-I design

Exceptional cutting performance

P-I Conical Drills' performance in dense bone, at the highest recommended rotation, without gradual diameter increments and applying constant feeding, present a very low friction therefore lower temperature transmission to bone tissue. Data on file.



Easy, simplified installation

Maximum of 3
low speed steps

(!) Except for MT-F Ø 4.8 Implant, 3 or 4 low speed steps are used. See Surgical Sequence. Spade and Round burr are optional. No pilot drill, counter sink or screw-tap required.



Insertion Driver


Handpiece • Manual • Torque Wrench





Biological Width

Driver reference



(!) The horizontal Implant Insertion Driver's mark is at approximately 3 mm and serves as a minimum Biological Width vertical reference for Implant platform positioning when completely covered by the lowest point of the soft tissue, the gingival margin. For further submersion, verify available prosthetic Component dimensions. Implant Insertion Driver dots, upper hexagon and flat areas are indexed to the Implant's hexagonal index and Sleeve • GS flat and slots. When used with Torque Wrench, the upper hexagonal portion of the Implant Insertion Driver should be entirely connected to the hexagon of the Driver Adapter.




 Guided &
conventional
SURGERY



Precision • Sensitivity • Versatility

Fully Guided Surgery novel stepwise approach



Drill frictionless
through sleeve

Potential reduction of both
temperature and debris
release on surgical site

Initial and Conical Drills do not rotate
against surgical guide or Sleeves




files for open
source software
piبرانemark.com



DentalCAD® exoplan





● MT-F
Implants

Platform Ø	3.3	3.5	3.9	4.6
h				
18		172319		
15	172297	172302	172384	
13	172296	172301	172383	172306
11.5	172295	172300	172382	172305
10	172294	172299	172381	172304
8.5	172293	172298	172380	172303
7		172318	172379	172321
6		172317	172378	172320



Implant Ø

3.3

3.75

4.1

4.8



Short
Implants

D DIGITAL
prosthetics




⦿ Components

Soft Tissue Healing Cover Screw • Healing Abutment



Biological Width

Concave or Parallel emergence Healing • Soft Tissue contouring
 Potential for more soft tissue volume
 Minimized cortical bone removal for sub-crestal Implants



One Cover Screw

For all Implants and Platforms • MT Interface



			N	R	W
			h		
Divergent		5.5	172545	172546	172547
		4.5	171199	171202	171205
		3	171198	171201	171204
		1.5	171197	171200	171203
Parallel			SELECTION		
		4.5	171190	171193	171196
		3	171189	171192	171195
		1.5	171188	171191	171194

Healing Abutment

Cover		171104
-------	---	--

Cover Screw

Regular Cover Screw is used.

Healing Abutment total height is +1.0mm higher than nominal margin height, except h = 5.5, + 1.5mm.



Conical Abutment

Single or multiple, screw-retained prosthesis



Biological Width

Concave emergence • Increased soft tissue volume

Minimizes cortical bone removal to install Abutment

Multiple and Single prosthesis

Universal Ø 4.8 Multi Unit platform

Straight Conical Abutment has double indexation

Select engaging components for single unit


DIGITAL

Provisional and CAD/CAM prosthesis

Adjustable post heights: 10, 9, 8, 7, 5.5 and 4.5mm



Direct Geometry for single and multiple unit • bridge



 **Conical Abutment**

- Regular Abutment is used.

- (!) Conical Abutment prosthetic Platform has $\varnothing$ 4.8 mm.

Cylinder • Conical Abutment ENG (Engaging) option is available for straight Conical Abutment only.

● Abutment Cemented Cylinder

Single or multiple, cement-retained prosthesis



Biological Width

Concave emergence • Potential for more soft tissue volume

Minimizes cortical bone removal to install Abutment

Anterior and posterior

Indexed • 6 and 4mm post heights

Single and multiple Castable Cemented Cylinders ○ ●




One-time one-abutment option

Prosthetic procedures over Abutment or Implant Platform

Zero margin Abutment

For limited interproximal spaces “0”



Cylinder



N

R

W

Non-Engaging

Castable 6mm (L)

161413

161418

161423

Castable 4mm

161463

101747

101977

Engaging

Castable 6mm (L)

161414

161419

161424

Castable 4mm

161464

101746

101976

Analog		6mm (L)	161410	161415	161420
		4mm	161462	101745	101975

Impression Coping		Closed Tray, 6mm (L)	161412	161417	161422
		Closed Tray, 4mm	161461	101744	101974

Healing Cap		6mm (L)	161411	161416	161421
		4mm	161460	101743	101973

4mm		h			
		4.5	171157	171162	171167
		3	171156	171161	171166
		1.5	171155	171160	171165
		0.8	171154	171159	171164

6mm • Long (L)		4.5	171142	171147	171152
		3	171141	171146	171151
		1.5	171140	171145	171150
		0.8	171139	171144	171149
		"0"	171138	171143	171148

Abutment Cemented Cylinder

(!) "0" margin Abutment is not compatible with Healing Cap, Impression Coping or Cylinder.



Contour & Esthetic Abutments

Single or multiple, cement-retained prosthesis



Increased Biological Width

Concave emergence • Potential for more soft tissue volume

Minimizes cortical bone removal to install Abutment

Robust design

Adjustable • Straight and 17°

Delicate slim profile

Adjustable • Straight and 15°

Impression at Implant Platform


Short and Long Impression Coping

Open and Closed Tray

Contour



Esthetic



			N	R	W
			h		
Contour 17°		4.5	171116	171119	171122
		3	171115	171118	171121
		1.5	171114	171117	171120
Contour Straight		4.5	171107	171110	171113
		3	171106	171109	171112
		1.5	171105	171108	171111

Contour Abutment

		h			
Esthetic 15°		4.5	171178	171181	●
		3	171177	171180	●
		1.5	171176	171179	●
Esthetic Straight		4.5	171171	171175	●
		3	171170	171174	●
		1.5	171169	171173	●
		0.8	171168	171172	●

Esthetic Abutment

Implant Analog		●	171212	●	
Impression Coping Implant		Open Tray	●	171206	●
		Closed Tray	●	171209	●
		Open Tray, Long	●	172418	●
		Closed Tray, Long	●	172417	●

Implant Impression



Cylinders over Implant

Single or multiple, cement or screw-retained prosthesis



ENHANCED

Ceramic Surface



DIGITAL
prosthetics

Provisional, conventional and CAD/CAM

Platform diameter and margin height options

Adjustable post heights: 6, 5 and 4mm



Increased Biological Width

Concave emergence • Potential for more soft tissue volume

Minimizes cortical bone removal to install Cylinder

Definitive • Overcasting

Main body CoCrMo and waxing sleeve POM with retentions

Provisional and Custom Healing • Titanium

Diameter options • Flat areas and trapezoidal retentions

Impression at Implant Platform

Short and Long, Open and Closed Tray, Impression Coping



			N	R	W	
Cylinder DIGITAL and Conventional		○ Non-Engaging	h			
			4.5	172534	172536	172541
			3	172533	172539	172544
			1.5	172532	172538	172543
		0.8	172531	172537	172542	
		● Engaging	4.5	172401	172405	172409
			3	172404	172408	172412
			1.5	172403	172407	172411
			0.8	172402	172406	172410
		Cylinder Conventional		○ Non-Engaging	Titanium	172530
● Engaging						
CoCrMo	171183			171185	171187	
Titanium	171182			171184	171186	

○ Cylinder over Implant

Implant Analog		●	171212	●	
Implant Impression Coping		Open Tray	●	171206	●
		Closed Tray	●	171209	●
		Open Tray, Long	●	172418	●
		Closed Tray Long	●	172417	●

Implant Impression

D

 See DIGITAL Prosthetics.
 ● Regular Components are used.
 (!) Cylinder over Implant • Ti with height options can also be used to support prosthesis manufactured by CAD/CAM.

D DIGITAL prosthetics



Manufacturing
Prosthetics • 3D Model



Components
and Direct Geometries



Design
Libraries



Scan
Intraoral • Desk



Scan Body
Implant • Conical Abutment



files for open
source software
pibranemark.com

3shape

exocad





Link CAD | CAM



h			
4.5	●	171137	●
3	●	171136	●
1.5	●	171135	●
0.8	●	171134	●

Link C



	<u>Small</u>	<u>Large</u>	
4.5	172422	172426	●
3	172421	172425	●
1.5	172420	172424	●
0.8	172419	172423	●

Link

Implant • MT
for cast or 3D models



●	172416	●
---	--------	---

Analog + Digital

Implant • MT



PEEK	●	161469	●
Titanium	●	161563	●

Scan Body

● MT interface is the same for all platforms. Scan Body and Analog + Digital are used.
(!) Links C 3.5 • Small and 4.1 • Large have similar post geometries to Dentsply Sirona, CEREC®.

D Cylinder over Implant

ENHANCED
Ceramic Surface

N R W

Cylinder
DIGITAL and
Conventional



h			
	Non-Engaging		
4.5	172534	172536	172541
3	172533	172539	172544
1.5	172532	172538	172543
0.8	172531	172537	172542
	Engaging		
4.5	172401	172405	172409
3	172404	172408	172412
1.5	172403	172407	172411
0.8	172402	172406	172410

Cylinder over Implant

Implant • MT
for cast or 3D models



• 172416 •

Analog + Digital

Implant • MT



PEEK	•	161469	•
Titanium	•	161563	•

Scan Body

D Cylinder • Conical Abutment

N

R

W

Titanium
DIGITAL and
Conventional



Non-Engaging

101142



Engaging

171248

Cylinder

Conical Abutment
for cast or 3D models



172453

Analog + Digital

Conical Abutment



Non-Engaging



PEEK

161471

Titanium

161565

Engaging



PEEK

161556

Titanium

161566

Scan Body

Milling Blank



A - Arum



172443



Milling Blank

Analog Fixation



Conical Abutment Tool

70170003

Implant MT Tool

70170002

Screw

70160502

Digital Analog
with fixation tool



Conical abutment, NENG

70150007

Implant MT

70150005

Scan Body
Titanium



Conical abutment, NENG

80620277

Implant MT Interface



80620276



Third Party Acessories

● MT interface is the same in all platforms. Scan Body, Analog + Digital and Milling Blank are used.


■ Manufactured by CoreMedTech. Check latest CoreMedTech libraries version and availability.

(!) Cylinder • Conical Abutment ENG (Engaging) option is available for straight Conical Abutment only.



D Dimensions and References

Link



Cylinder over Implant



Cylinder • Ti • Conical Abutment



Direct Geometry

Cylinder
Conical
Abutment



Engaging



Non-Engaging

Bridge
over
Implant



(!) Direct Geometries are used to manufacture direct screw-retained CAD/CAM prosthesis without P-I titanium base (Links or Cylinders), using spare P-I Screws. Always check latest library version and availability.



Kit

Conventional



width	254 mm
height	40
depth	130

Advanced Stainless Steel

181036

Kit

Tray options




w	202
h	67
d	158

Advanced
Polymer
181022



w	120
h	40
d	80

Advanced
Compact
181023



w	120
h	40
d	80

Prosthetic
Compact
181029

(!) Prosthetic Kit does not include Torque Wrench and Adapters.
(!) Conventional Kit Composition at pibranemark.com
(!) Reference number for ordering purposes only. Instruments and Tray delivered separately.

Instruments • Conventional

Drills				Ø	
	Initial			2.2	141138
				2.8	141146
	Conical			3.4	141148
				3.8	141314
				4.6	141152
				4.8	141315
	Dense			3.3	141213
				3.75	141316
				4.0	141215
				4.8 5.0	141317

Implant Insertion Driver				Medium	Long
	All Systems			131139	131140

Tools					
	Guide Pin			2.2 2.8	131114
				2.8 3.8	131115
				2.2 2.8 C	141535
				2.8 3.8 C	141536
	Drill Extension				131028
	Spade			1.5	141319 ▲
	Round Burr			1.8	141001 ▲

▲ Optional Instruments not included in the Advanced Kit contents.

 Guided &
conventional
SURGERY



Wash Tray • GS
Surgical + Prosthetic



Guided & Conventional

181040

All Instruments



Guided +

181041

Guided Surgery Instruments only
includes Prosthetic Instruments



Tray • GS

131162

Tray only

width	190
height	61
depth	138



Stop Organizer • GS

141399

Organizer [9] only

Instruments • Guided & Conventional

Drills




Crestal Drill

Narrow

Regular

141330

141335



Initial

Conical

			
Ø	N • Long	R • Short	R • Long
2.2	141378	141336	141340
2.8	141381	141337	141341
3.4		141338	141342
3.8		141339	141343

Stop • GS Offset [9]


Set • 11 Stops - GS 181043
without organizer

			
h	9/10	9/10	9/11.5
6		141365	
7		141366	
8.5		141367	
10	141360	141368	
11.5	141361		141369
13	141362		141370
15	141363		141371

Sleeve




141375




141376

Implant Insertion Driver




131153

131152



Insertion Driver must be completely connected into the Implant internal hexagonal indexation. When used manually or with Torque Wrench, the upper hexagonal portion of the Insertion Driver must be entirely connected to the hexagon of the Driver Adapter.


Guide Fixation

	Sleeve - Pin	141382
	Pin Drill Ø1.7	141377
	Fixation Pin	131167



GUIDE STABILITY
Conical seating of Pin

Surgical Accessories

		Narrow	Regular
Punch			
	Soft Tissue Punch	141331	141332
Dense Drill			
	Dense Drill	141344	141345


Torque Wrench Kit



Prosthetic Instruments

	Hexagonal Ø 1.2	Short	131010
		Medium	131011
		Long	131012
	Conical Abutment Ø 2.0	Short	131016
		Medium	131017
	Handling Tool, Angled		102964
	Retriever MT	Short	141564
		Medium	131131

Accessories



		Depth Probe • GS [9] Offset tip + Round tip	141439	●
		Stop Removal Tool	131172	●

■ Manufactured by Elos MedTech Pinol A/S . Torque Wrench Kit includes Surgical & Prosthetic Adapters.
● Stop Removal Tool should be used when Surgical Adapter is not sufficient for the removal of Stop from Drill.
Depth Probe is not part of Kit • GS and should be ordered separately.

Stop • Conventional

		 R • Short	 R • Long
h			
	6	141592	
	7	141593	
	8.5	141594	
	10	141590	
	11.5	141591	141584
	13	141600	141586
	15	141601	141586
	18		141602

Trays



Polymer
131134



Stainless Steel
131117




Prosthetic
131135

Compact
131138

Lab & Replacement
Conical Abutment



Screw Cylinder
101984



Screw Impression Coping
Open Tray 101737
Closed Tray 171260

Implant



Coated Screw MT
171239



Screw Impression Coping
Open Tray 171207

Resonance Frequency Analysis

Penguin^{RFA} ⁽¹⁾



ISQ

Kit
55002-S ▲



MulTipeg

MT-F • Type 38

55046

Conical Abutment • Type 72

55080

Locator® Overdenture

Abutment ⁽²⁾



	N	R	W
h			
4	●	2203	●
3	●	2202	●
2	●	2201	●
1	●	2200	●

 (1) Manufactured by Integration Diagnostics Sweden AB .

▲ Penguin RFA Kit includes instrument, charger, MulTipeg driver and user's manual.

 (2) Manufactured by Zest Dental . Components and instruments not included in the P-I Catalog.

● MT interface is the same in all platforms. Regular Abutments are used.


P-I pegs are available from Osstell. MT-F Implant SmartPeg, type 21.

● MT-F

SURGICAL SEQUENCE



	rpm	600 - 1,200 Lowest possible rpm
	In-Out	△ Coordinated in-and-out movement of Drills for better cooling
	ITV	≤ 70 Ncm Insertion Torque Value
	Full Length	Prepare at planned full length of Implant position combining the Drills GS and Stops GS in accordance to planning
	Irrigation	Constant irrigation to the insertion margin of Conical Drills



△ Drills are less than 1 mm longer than Drill marks

SURGICAL SEQUENCE

 3.3	<table><tr><th>2.2</th><th>2.8</th><th>3.4</th></tr><tr><td>✓</td><td>S ↓ N</td><td>D ↓</td></tr></table>	2.2	2.8	3.4	✓	S ↓ N	D ↓	 3.75	<table><tr><th>2.2</th><th>2.8</th><th>3.4</th></tr><tr><td>✓</td><td>S</td><td>N ↓ D</td></tr></table>	2.2	2.8	3.4	✓	S	N ↓ D				
2.2	2.8	3.4																	
✓	S ↓ N	D ↓																	
2.2	2.8	3.4																	
✓	S	N ↓ D																	
 4.1	<table><tr><th>2.2</th><th>2.8</th><th>3.4</th><th>3.8</th></tr><tr><td>✓</td><td>S</td><td>N Optional for Dense</td><td>D</td></tr></table>	2.2	2.8	3.4	3.8	✓	S	N Optional for Dense	D	 4.8	<table><tr><th>2.2</th><th>2.8</th><th>3.8</th><th>4.6</th></tr><tr><td>✓</td><td>✓</td><td>S</td><td>N ↓ D</td></tr></table>	2.2	2.8	3.8	4.6	✓	✓	S	N ↓ D
2.2	2.8	3.4	3.8																
✓	S	N Optional for Dense	D																
2.2	2.8	3.8	4.6																
✓	✓	S	N ↓ D																

Bone Density

S Soft N Normal D Dense


 6 mm
Drill tip

(!) The subsequent Conical Drill, in terms of diameter, should be considered with a drilling depth of 6 mm, in order to not exceed 70 Ncm of insertion torque value. The use of Dense Drills (15 – 50 rpm) can also be considered to lower the insertion torque value. Soft Tissue Punch (15 – 50 rpm).


🔗 See Guided Surgery Surgical Sequence at pibranemark.com. Guided Surgery Drills consider a [9mm] offset and, when used with Drill stops, allow for limiting the total length of osteotomy with the objective of providing predetermined Drill length and orientation through the surgical guide. Height repositioning for Sleeve and Stop selection required for [10.5] Offset.

One interface

for all Implant dimensions



Dimensions



(!) Images are for illustrative purposes only. Measurements in millimeters.

One Driver

Manual • Torque Wrench



 1.2



Torques

Ncm

MT-F Implants

≤ 70

Abutments

Cylinders over Implant

25

Links, Milling Blank

Cylinders • Conical Abutment

15

Cover Screw


Healing Abutments

Impression Copings

Manual

Scan Bodies

(!) Recommended Torques. Abutment and Components torques should not exceed the torque obtained at Implant installation.
(!) One Prosthetic Driver, except Locator® and straight Conical Abutment.



LifeTime GUARANTEE





Developed By P-I Brånemark

MT-F SMART GUIDE-009-ENG-20.11.2024
P-I Brånemark. All rights reserved.

pibranemark.com

 SIC invent Sweden AB
Stora Åvägen 21
43634 Askim, Sweden

CE 2460 CE