



Guided &
conventional
SURGERY



Developed By P-I Brånemark




DIGITAL



Fully Guided Surgery

novel stepwise approach



Precision • Sensitivity • Versatility



Guided &
conventional



DIGITAL
prosthetics

The P-I Brånemark • Guided & Conventional novel stepwise approach was developed to simplify clinical procedures and reduce the potential risks associated with guided surgery providing flexibility for conventional procedures by using the same Instruments and Kit.

Designed to increase tridimensional accuracy of Implant positioning, the Instruments overcome the inherent tactile feedback reduction, visualization and access limitations of surgical site preparation and Implant installation when using surgical guides.

Combined with the MT-F Implant System, the novel Guided & Conventional concept is a stepwise approach to the modern, immediate and digitalized implantology which keeps clinicians focused on patients and in control of procedures, even when plans change during surgery.

The P-I Guided & Conventional provides a predictable, minimally invasive and versatile solution for precise tridimensional positioning of Implants with an optimal esthetic outcome when used with integrated digital diagnostic, planning and prosthetic software.



Contents

Guided & Conventional

Novel Stepwise Approach Technology	5
Combinations	10
Dimensions	11
Alignments	13

Kit • Guided & Conventional Surgery	14
--	----

Instruments	15
--------------------	----

MT-F Guided Surgical Sequence	19
--------------------------------------	----





 Guided & conventional

Wash Tray
Surgical + Prosthetic





State-of-the-art

Less friction

Only titanium Stop is in contact with Sleeve

Less trauma

Constant apical angle • Diamond Like Carbon

High performance

P-I custom technology




Exceptional cutting performance

P-I Conical Drills' performance in dense bone, at the highest recommended rotation, without gradual diameter increments and applying constant feeding, present a very low friction coefficient therefore lower temperature transmission to bone tissue. Data on file.



Drill frictionless through Sleeve


🔄 Screw-retained Stops are the only contact with Sleeves • Titanium to titanium
Initial and Conical Drills do not rotate against surgical guide or Sleeves
Potential reduction of both temperature and debris release on surgical site



🔄 Fully Guided • Precision in all steps

Crestal Drills are continuously guided through Sleeves in Offset [9] and [10.5]
Initial and Conical Drills are oriented by previous drilling steps and Stops

🔄 Screw-retained Stops provide accuracy and stability



Same Instruments for Guided & Conventional

Guided Surgery - GS identification

For conventional procedure

-  Guided • All
-  Narrow - N
-  Regular - R • Long - L
-  Regular - R • Short - S



Short and Long
Drills GS



Insertion
Drivers GS

Unmatched versatility


Less Instruments • Keyless • Mountless

Shorter Drills and Stops can be used for incremental preparation depth of osteotomy

Stops can be pre-assembled according to planning


▼ Incremental depth option

MT-F 3.75 x 10 • Offset [9]




Full depth option

MT-F 3.75 x 10 • Offset [9]



Surgical Sequence according to bone density

Implant installation



Effective Insertion Torque Value ITV ↗

Reduced contact area • Diamond Like Carbon • Less friction on Sleeve
Safe pick-up, installation and removal • Avoids undesired movement
Possible further Implant submersion through Sleeve

Handpiece • Manual • Torque Wrench

Combinations


Narrow


Regular
Short


Regular
Long


Drill • GS



Stop • GS 



Sleeve • GS



MT-F

[10.5] Offset [10.5] and Biological Width oriented positioning should consider drilling depth compensation. Offset [10.5] cannot be used with MT-F h = 15.



h=18 For MT-F h=18 use Drill Stop - GS - 15, and the final drilling depth is determined conventionally without surgical guide.

4.8 For MT-F Ø4.8 use Conical Drill - GS - 3.8 for soft bone. The following drilling steps are performed conventionally without surgical guide. Implant insertion can be performed through Sleeve - GS - R.



Consult Instructions for Use and Guided Surgical Sequence.
Measurements in millimeters. Images for illustrative purposes only.
Verify available libraries and updates before planning.



Dimensions


Drill • GS



Stop • GS



Sleeve • GS



d	1.8	3.5	5
D	4.2	4	6
h	0.7	1	1
H	7	4	4
f	5	4.5	6
F	5	4.5	7




Sleeve - GS - N can also be used in limited spaces and for marking the initial drilling steps preceding conventional (free hand) surgical site finalization for the sleeveless installation of MT-F Ø3.75, 4.1 and 4.8.



Consult Instructions for Use, Implant Dimensions and Guided Surgical Sequence.
Measurements in millimeters. Images for illustrative purposes only.
Verify available libraries and updates before planning.




MT-F Dimensions



One Interface for all Implants

D DIGITAL
prosthetics

Alignments




Flat to Buccal



- The flat part of the Sleeve - R and the Sleeves' slots can be positioned buccally to align with the prosthetic orientation of the Implant Insertion Driver (flat areas, upper hexagon and precision dot) which are aligned to Implant and Prosthetic Components' indexation.

Corner to Buccal





- For angled Conical Abutments, and in case of mesiodistal space limitation, the flat part of the Sleeve - R can be mesially oriented leading to a Sleeve corner to buccal positioning.
- When the Sleeve - R is positioned corner to buccal, the prosthetic orientation of the Implant Insertion Driver can be directed to the corner of the Sleeve - R, and between the Sleeves' slots, for a buccal indexation of Implant and Prosthetic Components (except with angled Conical Abutments).



⚠ A minimum distance of 3mm should be ensured between 2 Implants

Depth verification

Offset [9] example



Wash Tray • GS

Surgical + Prosthetic



Guided & Conventional

181040

All Instruments



Guided +

181041

Guided Surgery Instruments only
includes Prosthetic Instruments



Tray • GS

131162

Tray only

width	190
height	61
depth	138



Stop Organizer • GS

141399

Organizer [9] only

Instruments • GS

Drill




Crestal Drill

Narrow

Regular

141330

141335



Initial

Conical

			
Ø	N • Long	R • Short	R • Long
2.2	141378	141336	141340
2.8	141381	141337	141341
3.4		141338	141342
3.8		141339	141343

Stop • GS Offset [9]

Set • 11 Stops - GS 181043
without organizer

h			
6		141365	
7		141366	
8.5		141367	
10	141360	141368	
11.5	141361		141369
13	141362		141370
15	141363		141371

Sleeve




141375




141376

Implant Insertion Driver



131153

131152



Insertion Driver must be completely connected into the Implant internal hexagonal indexation. When used manually or with Torque Wrench, the upper hexagonal portion of the Insertion Driver must be entirely connected to the hexagon of the Driver Adapter.




Guide Fixation

	Sleeve - Pin	141382
	Pin Drill Ø1.7	141377
	Fixation Pin	131167



GUIDE STABILITY
Conical seating of Pin

Surgical Accessories

		Narrow	Regular
Punch			
	Soft Tissue Punch	141331	141332
Dense Drill			
	Dense Drill	141344	141345

Torque Wrench Kit




C14569

HIGH TORQUE VERSION
originally developed for OSPOL • P-I

Prosthetic Instruments

	Hexagonal Ø 1.2	Short	131010
		Medium	131011
		Long	131012
	Conical Abutment Ø 2.0	Short	131016
		Medium	131017
	Handling Tool, Angled		102964
	Retriever MT	Short	141564
		Medium	131131

Surgical Accessories • GS



In the event Stop cannot be easily removed by hand, use handpiece in reverse or the Surgical Adapter, with or without Torque Wrench, to stabilize the Drill and apply manual counter torque on the Stop using the Stop Removal Tool.

	Depth Probe • GS [9] Offset tip + Round tip	141439	●
	Stop Removal Tool	131172	●


■ Manufactured by Elos MedTech Pinol A/S . Torque Wrench Kit includes Surgical & Prosthetic Adapters.
● Stop Removal Tool should be used when Surgical Adapter is not sufficient for the removal of Stop from Drill.
Depth Probe is not part of Kit • GS and should be ordered separately.

Stop • Conventional



STOP
conventional
FREE
HAND

Short • Long
Drills GS



R • Short

R • Long

h		
6	141592	
7	141593	
8.5	141594	
10	141590	
11.5	141591	141584
13	141600	141586
15	141601	141586
18		141602

MT-F Guided Surgical Sequence



rpm

600 - 1,200
Lowest possible rpm



In-Out

△ Coordinated in-and-out movement of Drills for better cooling



ITV

≤ 70 Ncm
Insertion Torque Value



Full Length

Prepare at planned full length of Implant position combining the Drills GS and Stops GS



Irrigation

Constant irrigation to the insertion margin of Conical Drills

△ Drills are less than 1 mm longer than Drill marks



3.3

Crestal	2.2	2.8	3.4
✓	✓	S ↓	D ↓
		N	

✓ ✓ S ↓ D ↓
N



3.75

Crestal	2.2	2.8	3.4
✓	✓	S	N ↓
			D

✓ ✓ S N ↓
D



4.1

Crestal	2.2	2.8	3.4	3.8
✓	✓	S	N	D
		 Optional for Dense		

✓ ✓ S N D
..... Optional for Dense



4.8

Crestal	2.2	2.8	3.8	4.6
✓	✓	✓	S	N ↓
				D

✓ ✓ ✓ S N ↓
D

▼ Incremental depth option

Bone Density

S Soft **N** Normal **D** Dense

↓ 6 mm Drill tip



The subsequent Conical Drill, in terms of diameter, should be considered with a drilling depth of 6 mm, in order to not exceed 70 Ncm of Insertion Torque Value. The use of Dense Drills (15 - 50 rpm) can also be considered to lower the Insertion Torque Value. The use of Soft Tissue Punch (15 - 50 rpm) precedes Crestal Drill - GS use and can also be used in conventional surgery procedures.



Consult Instructions for Use and Guided Surgical Sequence in detail.





Guided Surgical Sequence

in detail

MT-F		Guide		Guided Surgical Sequence						
Ø	h	Offset	Sleeve	Punch optional 15-50 rpm	1 Crestal 600 - 1,200 rpm	2 Initial 600 - 1,200 rpm	3 Bone Density 600 - 1,200 rpm			Implant Driver ≤ 70 Ncm
							SOFT	NORMAL	DENSE	
3.3	8.5	[9]	R	R	R	Ø 2.2 + 8.5	Ø 2.8 + 6 ▲	Ø 2.8 + 8.5	Ø 3.4 + 6	R
3.3	10.0	[9]	R	R	R	Ø 2.2 + 10	Ø 2.8 + 6 ▲	Ø 2.8 + 10	Ø 3.4 + 6	R
3.3	11.5	[9]	R	R	R	Ø 2.2 + 11.5	Ø 2.8 + 6 ▲	Ø 2.8 + 11.5	Ø 3.4 + 6	R
3.3	13.0	[9]	R	R	R	Ø 2.2 + 13	Ø 2.8 + 6 ▲	Ø 2.8 + 13	Ø 3.4 + 6	R
3.3	15.0	[9]	R	R	R	Ø 2.2 + 15	Ø 2.8 + 6 ▲	Ø 2.8 + 15	Ø 3.4 + 6	R
3.75	6.0	[9]	R	R	R	Ø 2.2 + 6	Ø 2.8 + 6	Ø 3.4 + 6 ▲	Ø 3.4 + 6	R
3.75	7.0	[9]	R	R	R	Ø 2.2 + 7	Ø 2.8 + 7	Ø 3.4 + 6 ▲	Ø 3.4 + 7	R
3.75	8.5	[9]	R	R	R	Ø 2.2 + 8.5	Ø 2.8 + 8.5	Ø 3.4 + 6 ▲	Ø 3.4 + 8.5	R
3.75	10.0	[9]	R	R	R	Ø 2.2 + 10	Ø 2.8 + 10	Ø 3.4 + 6 ▲	Ø 3.4 + 10	R
3.75	11.5	[9]	R	R	R	Ø 2.2 + 11.5	Ø 2.8 + 11.5	Ø 3.4 + 6 ▲	Ø 3.4 + 11.5	R
3.75	13.0	[9]	R	R	R	Ø 2.2 + 13	Ø 2.8 + 13	Ø 3.4 + 6 ▲	Ø 3.4 + 13	R
3.75	15.0	[9]	R	R	R	Ø 2.2 + 15	Ø 2.8 + 15	Ø 3.4 + 6 ▲	Ø 3.4 + 15	R
3.75	18.0	[9]	R	R	R	Ø 2.2 + 15	Ø 2.8 18 mark	Ø 3.4 6 mark ▲	Ø 3.4 18 mark	R
4.1	6.0	[9]	R	R	R	Ø 2.2 + 6	Ø 2.8 + 6	Ø 3.4 + 6 ▲	Ø 3.8 + 6	R
4.1	7.0	[9]	R	R	R	Ø 2.2 + 7	Ø 2.8 + 7	Ø 3.4 + 7 ▲	Ø 3.8 + 7	R
4.1	8.5	[9]	R	R	R	Ø 2.2 + 8.5	Ø 2.8 + 8.5	Ø 3.4 + 8.5 ▲	Ø 3.8 + 8.5	R
4.1	10.0	[9]	R	R	R	Ø 2.2 + 10	Ø 2.8 + 10	Ø 3.4 + 10 ▲	Ø 3.8 + 10	R
4.1	11.5	[9]	R	R	R	Ø 2.2 + 11.5	Ø 2.8 + 11.5	Ø 3.4 + 11.5 ▲	Ø 3.8 + 11.5	R
4.1	13.0	[9]	R	R	R	Ø 2.2 + 13	Ø 2.8 + 13	Ø 3.4 + 13 ▲	Ø 3.8 + 13	R
4.1	15.0	[9]	R	R	R	Ø 2.2 + 15	Ø 2.8 + 15	Ø 3.4 + 15 ▲	Ø 3.8 + 15	R

▼ Incremental depth option

Possible use of shorter Stops - GS preceding final depth for increased angular accuracy through Sleeve

R = Regular

● Short ● Long

▲ Optional step for dense bone



Guided Surgical Sequence

in detail

MT-F		Guide		Guided Surgical Sequence						
Ø	h	Offset	Sleeve	Punch optional	1 Crestal	2 Initial	3 Bone Density			Implant Driver
				15-50 rpm	600 - 1,200 rpm	600 - 1,200 rpm	SOFT	NORMAL	DENSE	≤ 70 Ncm
4.8	6.0	[9]	R	R	R	Ø 2.2 + 6	Ø 2.8 + 3.8 + 6	Ø 4.6 6 mark ▲	Ø 4.6 6 mark	R
4.8	7.0	[9]	R	R	R	Ø 2.2 + 7	Ø 2.8 + 3.8 + 7	Ø 4.6 6 mark ▲	Ø 4.6 7 mark	R
4.8	8.5	[9]	R	R	R	Ø 2.2 + 8.5	Ø 2.8 + 3.8 + 8.5	Ø 4.6 6 mark ▲	Ø 4.6 8.5 mark	R
4.8	10.0	[9]	R	R	R	Ø 2.2 + 10	Ø 2.8 + 3.8 + 10	Ø 4.6 6 mark ▲	Ø 4.6 10 mark	R
4.8	11.5	[9]	R	R	R	Ø 2.2 + 10	Ø 2.8 + 3.8 11.5 mark	Ø 4.6 6 mark ▲	Ø 4.6 11.5 mark	R
4.8	13.0	[9]	R	R	R	Ø 2.2 + 10	Ø 2.8 + 3.8 13 mark	Ø 4.6 6 mark ▲	Ø 4.6 13 mark	R
3.3	10.0	[9]		N	N	Ø 2.2 + 10	Dense N ▲	Ø 2.8 + 10	Dense N	N
3.3	11.5	[9]	N	N	N	Ø 2.2 + 11.5	Dense N ▲	Ø 2.8 + 11.5	Dense N	N
3.3	13.0	[9]	N	N	N	Ø 2.2 + 13	Dense N ▲	Ø 2.8 + 13	Dense N	N
3.3	15.0	[9]	N	N	N	Ø 2.2 + 15	Dense N ▲	Ø 2.9 + 15	Dense N	N

Incremental depth option

Possible use of shorter Stops - GS preceding final depth for increased angular accuracy through Sleeve

R = Regular N = Narrow

• Short • Long • Narrow

▲ Optional step for dense bone



Full

Prepare at planned full length of Implant position combining the Drills - GS and Stops - GS



In-Out

△ Coordinated in-and-out movement of Drill for better cooling



Irrigation

Constant cooled irrigation to the insertion margin of Drills



6 mm Drill tip

The subsequent Conical Drill, in terms of diameter, should be considered with a drilling depth of 6 mm, in order to not exceed 70 Ncm of Insertion Torque Value. The use of Dense Drills (15 - 50 rpm) can also be considered to lower the Insertion Torque Value.



[10.5]

Offset [10.5] and Biological Width oriented positioning should consider drilling depth compensation



Consult Instructions for Use



MT-F Guided Surgical Sequence in detail






Developed By P-I Brånemark

P-I GUIDED & CONVENTIONAL-001-ENG-1710-2024
P-I Brånemark. All rights reserved.

pibranemark.com

 SIC invent Sweden AB
Stora Åvägen 21
43634 Askim, Sweden



files for open
source software
pibranemark.com

exocad DentalCAD® exoplan

